

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge CB4 0WY, United Kingdom Email: info@biorbyt.com

Phone: +44 (0)1223 859 353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

369 Pine Street, Suite #103, San Francisco California, 94104, United States Email: info@biorbyt.com

Phone: +1 (415) 906 5211 | Fax: +1 (415) 651 8558

Product Datasheet

Dog ADM(Adrenomedullin) ELISA Kit orb1994871

Description Dog ADM(Adrenomedullin) ELISA Kit

Reactivity Canine

Range 15.63-1000 pg/mL

Sensitivity 4.42 pg/mL

Note For research use only.

Sample Types Serum, plasma and other biological fluids.

Expiration Date Please enquire.

Application Notes The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate

provided in this kit has been pre-coated with an antibody specific to Dog ADM. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog ADM. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog ADM, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ± 10nm. The concentration of Dog ADM in the samples is then

determined by comparing the OD of the samples to the standard curve.