

Product Datasheet

MYOGLOBIN Method: latex agglutination orb668846

Description	Polystyrene latex beads coated with antibodies to human myoglobin in the presence of myoglobin in human serum form visible agglutination. The rate of agglutination correlates with myoglobin concentration in serum
Hazard Information	Non-Toxic
Target	MYOGLOBIN Method: latex agglutination
Storage	Store at 2-8°C
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Increased levels associated with heart diseases. Myoglobin is the first marker of infarction.