

Product Datasheet

PROTEIN (TOTAL) Method: colorimetric, endpoint (biuret) orb668847

Description	The copper ions in alkaline biuret reagent react with peptide bonds of serum proteins to form a purple color. The intensity of the color is proportional to the total protein concentration
Hazard Information	Non-Toxic
Target	PROTEIN (TOTAL) Method: colorimetric, endpoint (biuret)
Storage	Store at 2-8°C
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Increased levels observed in cases of dehydration.