

Biorbyt Ltd.

5 Orwell Furlong, Cowley Road, Cambridge CB4 0WY, United Kingdom Email: info@biorbyt.com

Phone: +44 (0)1223 859 353 | Fax: +44 (0)1223 280 240

Biorbyt LLC.

369 Pine Street, Suite #103, San Francisco California, 94104, United States Email: info@biorbyt.com

Phone: +1 (415) 906 5211 | Fax: +1 (415) 651 8558

## **Product Datasheet**

## PROTEIN (strips for protein detection) orb668849

**Description** Principle:Albumin in urine reacts with dye on strip bring to color change.

Color intensity is correlated with amount of protein. The strips are applied

for semi quantitative determination of protein in urine

Hazard Information Non-Toxic

Target PROTEIN (strips for protein detection)

Storage Room Temperature

**Note** For research use only.

**Expiration Date** 12 months from date of receipt.

**Application Notes** Protein in urine is indicator of kidney dysfunction