

## Product Datasheet

### UREA Method : colorimetric, endpoint orb668852

Description	Urea + H <sub>2</sub> O+urease= NH <sub>3</sub> + CO <sub>2</sub> nitroprusside NH <sub>3</sub> + salicilate-hypochlorite = colored complex
Hazard Information	Non-Toxic
Target	UREA Method : colorimetric, endpoint
Storage	Store at 2-8°C
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Increased levels associated with renal diseases as well as dehydration, diabetic coma.