

Product Datasheet

URIC ACID Method: ultraviolet, endpoint orb668853

Description	Uric acid + H ₂ O + O ₂ + uricase = allantoin + CO ₂ + H ₂ O ₂ Absorbance at 293 nm is measured before and after reaction. Differences in absorbances show concentration of uric acid
Hazard Information	Non-Toxic
Target	URIC ACID Method: ultraviolet, endpoint
Storage	Store at 2-8°C
Note	For research use only.
Expiration Date	12 months from date of receipt.
Application Notes	Increased levels associated with gout, leukemia, severe renal impairment.