

Product Datasheet

Mouse Anti-Rabbit IgG antibody orb973625

Description	Mouse Anti-Rabbit IgG antibody
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	Native Rabbit IgG
Form/Appearance	Lyophilized powder
Preservatives	0.01M PBS (pH7.4)
Storage	Store at -20°C stable for 2 years (lyophilized powder). Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
Clone ID	4C11
Dilution Range	WB=1:5000-50000,IHC-P=1:100-500,IHC-F=1:1000-5000,ICC=1:1000-5000,IF=1:1000-5000,ELISA=1:5000-50000
Expiration Date	12 months from date of receipt.