

## Product Datasheet

### LDL antibody {FITC} orb6310-FITC

Description	Rabbit polyclonal antibody against LDL conjugated to FITC.
Species/Host	Rabbit
Reactivity	Human, Mouse
Predicted Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	Lipoprotein, low density from human plasma
Target	LDL {FITC}
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	N/AkDa
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.