

Product Datasheet

IKK γ (Phospho-T23) antibody orb1151171

Description	Rabbit polyclonal antibody to IKK γ (Phospho-T23)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding T23 of human IKK γ .
Storage	Store at -20 $^{\circ}$ C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O15111
Entrez	1147
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.