

Product Datasheet

NF-kB p65/RelA antibody orb1150744

Description	Rabbit polyclonal antibody to NF-kB p65/RelA
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human NF-kB p65/RelA(NP_068810.3).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	65KDa
Uniprot ID	Q04206
NCBI	5970
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:100IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.