

Product Datasheet

LRP6 (Phospho-S1490) antibody orb1150137

Description	Rabbit polyclonal antibody to LRP6 (Phospho-S1490)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S1490 of human LRP6.
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	180KDa
Uniprot ID	O75581
Entrez	4040
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.