

Product Datasheet

Mu Opioid Receptor(MOR) antibody orb1150050

Description	Rabbit polyclonal antibody to Mu Opioid Receptor(MOR)
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Mu Opioid Receptor(MOR)(NP_000905.3).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	10-20kDa/33-55kDa
Uniprot ID	P35372
NCBI	4988
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.