

Product Datasheet

CXXC4 antibody orb1150048

Description	Rabbit polyclonal antibody to CXXC4
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CXXC4 (NP_079488.2).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9H2H0
Entrez	80319
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.