

Product Datasheet

LPHN2 antibody orb1150006

Description	Rabbit polyclonal antibody to LPHN2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human LPHN2 (NP_036434.1).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O95490
Entrez	23266
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.