

Product Datasheet

HRH1 antibody orb1149999

Description	Rabbit polyclonal antibody to HRH1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant protein of human HRH1.
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	60KDa
Uniprot ID	P35367
Entrez	3269
Dilution Range	WB 1:1000 - 1:2000IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.