

Product Datasheet

Beclin 1 antibody orb1149963

Description	Rabbit polyclonal antibody to Beclin 1
Reactivity	Human, Mouse
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human Beclin 1(NP_003757.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	52KDa
Uniprot ID	Q14457
NCBI	8678
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.