

Product Datasheet

PTPN22 antibody orb1149935

Description	Rabbit polyclonal antibody to PTPN22
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	Recombinant fusion protein of human PTPN22(NP_036543.4).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	20kDa/75kDa/78kDa/85kDa/88kDa/91kDa
Uniprot ID	Q9Y2R2
NCBI	26191
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.