

Product Datasheet

SHOX antibody orb1149929

Description	Rabbit polyclonal antibody to SHOX
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	Recombinant fusion protein of human SHOX(NP_000442.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	32kDa
Uniprot ID	O15266
NCBI	6473
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.