

Product Datasheet

CELSR1 antibody orb1149901

Description	Rabbit polyclonal antibody to CELSR1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human CELSR1 (NP_055061.1).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9NYQ6
Entrez	9620
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.