

Product Datasheet

ULK1 (Phospho-S638) antibody orb1149900

Description	Rabbit polyclonal antibody to ULK1 (Phospho-S638)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S638 of human Phospho-ULK1-S638.
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	150KDa
Uniprot ID	O75385
Entrez	8408
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.