

Product Datasheet

F2RL2 antibody orb1149892

Description	Rabbit polyclonal antibody to F2RL2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human F2RL2 (NP_004092.1).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	O00254
Entrez	2151
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.