

Product Datasheet

DKK3 antibody orb1149879

Description	Rabbit polyclonal antibody to DKK3
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human DKK3(NP_056965.3).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	60KDa
Uniprot ID	Q9UBP4
NCBI	27122
Dilution Range	WB 1:500 - 1:2000IF 1:50 - 1:100
Expiration Date	12 months from date of receipt.