

Product Datasheet

CSF1R (Phospho-Y809) antibody orb1149874

Description	Rabbit polyclonal antibody to CSF1R (Phospho-Y809)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y809 of human CSF1R (NP_005202.2).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P07333
Entrez	1436
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.