

Product Datasheet

CFC1 antibody orb1149835

Description	Rabbit polyclonal antibody to CFC1
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of Human CFC1.
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P0CG37
Entrez	55997
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.