

Product Datasheet

Parkin antibody orb1149832

Description	Rabbit polyclonal antibody to Parkin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human Parkin
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	52KDa
Uniprot ID	O60260
Entrez	5071
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200
Expiration Date	12 months from date of receipt.