

Product Datasheet

Methyl-Methyl-RhoA-H105 antibody orb1149827

Description	Rabbit polyclonal antibody to Methyl-Methyl-RhoA-H105
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human Methyl-Methyl-RhoA-H105, 3-methyl (tau-methyl).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P61586
Entrez	387
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.