

Product Datasheet

EDA2R antibody orb1149822

Description	Rabbit polyclonal antibody to EDA2R
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human EDA2R (NP_068555.1).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9HAV5
Entrez	60401
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.