

Product Datasheet

PKC delta (Phospho-S645) antibody orb1149820

Description	Rabbit polyclonal antibody to PKC delta (Phospho-S645)
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around S645 of human PKC delta(NP_006245.2).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
MW	80kDa
Uniprot ID	Q05655
NCBI	5580
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.