

Product Datasheet

Cytokeratin 8 (Phospho-S74) antibody orb1149812

Description	Rabbit polyclonal antibody to Cytokeratin 8 (Phospho-S74)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding S74 of human Cytokeratin 8
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P05787
Entrez	3856
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.