

Product Datasheet

SARS-CoV-2 NSP9 antibody orb1149806

Description	Rabbit polyclonal antibody to SARS-CoV-2 NSP9
Reactivity	Virus
Tested Applications	WB
Immunogen	Recombinant fusion protein of coronavirus NSP9(YP_009742616.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	15-20KDa
Uniprot ID	P0DTD1
NCBI	43740578
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.