

Product Datasheet

ARHGEF2 antibody orb1149803

Description	Rabbit polyclonal antibody to ARHGEF2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human ARHGEF2 (NP_004714.2).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	120KDa
Uniprot ID	Q92974
Entrez	9181
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.