

Product Datasheet

CtBP1 antibody orb1149799

Description	Rabbit polyclonal antibody to CtBP1
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CtBP1(NP_001319.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	47kDa
Uniprot ID	Q13363
NCBI	1487
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.