

Product Datasheet

IL19 antibody orb1149785

Description	Rabbit polyclonal antibody to IL19
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human IL19.
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	Q9UHD0
Entrez	29949
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.