

Product Datasheet

LAT (Phospho-Y191) antibody orb1149778

Description	Rabbit polyclonal antibody to LAT (Phospho-Y191)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A phospho specific peptide corresponding to residues surrounding Y191 of human LAT
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	36KDa
Uniprot ID	O43561
Entrez	27040
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.