

## Product Datasheet

### [KO Validated] HP1 alpha/CBX5 antibody orb1149777

|                            |                                                                                |
|----------------------------|--------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to [KO Validated] HP1 alpha/CBX5                    |
| <b>Species/Host</b>        | Rabbit                                                                         |
| <b>Reactivity</b>          | Human, Rat                                                                     |
| <b>Conjugation</b>         | Unconjugated                                                                   |
| <b>Tested Applications</b> | IF, WB                                                                         |
| <b>Immunogen</b>           | Recombinant fusion protein of human HP1 alpha/HP1 alpha/CBX5 (NP_001120794.1). |
| <b>Storage</b>             | Store at -20?. Avoid freeze / thaw cycles.                                     |
| <b>Note</b>                | For research use only.                                                         |
| <b>Isotype</b>             | IgG                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                     |
| <b>Purity</b>              | Affinity purification                                                          |
| <b>MW</b>                  | 22kDa                                                                          |
| <b>Uniprot ID</b>          | <a href="#">P45973</a>                                                         |
| <b>Entrez</b>              | <a href="#">23468</a>                                                          |
| <b>Dilution Range</b>      | WB 1:500 - 1:2000IF 1:50 - 1:200                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                |