

Product Datasheet

LEF1 antibody orb1149776

Description	Rabbit polyclonal antibody to LEF1
Reactivity	Human, Mouse, Rat
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human LEF1(NP_057353.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	25-58kDa
Uniprot ID	Q9UJU2
NCBI	51176
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.