

Product Datasheet

EVA1A antibody orb1149763

Description	Rabbit polyclonal antibody to EVA1A
Reactivity	Human, Mouse, Rat
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human EVA1A(NP_115557.1).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	17KDa
Uniprot ID	Q9H8M9
NCBI	84141
Dilution Range	WB 1:500 - 1:2000IHC 1:50 - 1:200
Expiration Date	12 months from date of receipt.