

Product Datasheet

SHC1 (Phospho-Y427) antibody orb1149762

Description	Rabbit polyclonal antibody to SHC1 (Phospho-Y427)
Reactivity	Human, Mouse
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y427 of human SHC1(NP_892113.4).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
MW	50KDa/55KDa
Uniprot ID	P29353
NCBI	6464
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.