

Product Datasheet

IL12A antibody orb1149756

Description	Rabbit polyclonal antibody to IL12A
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	A synthetic peptide of human IL12A(NP_000873.2).
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	35KDa/75KDa
Uniprot ID	P29459
NCBI	3592
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.