

Product Datasheet

MUC1 (Phospho-Y1203) antibody orb1149747

Description	Rabbit polyclonal antibody to MUC1 (Phospho-Y1203)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic phosphorylated peptide around Y1203 of human MUC1 (P15941).
Storage	Store at -20?. Avoid freeze / thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
Uniprot ID	P15941
Entrez	4582
Dilution Range	WB 1:500 - 1:2000
Expiration Date	12 months from date of receipt.