

## Product Datasheet

### Thrombin light chain antibody {HRP} orb158601-HRP

Description	Rabbit polyclonal antibody against Thrombin light chain conjugated to HRP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine,Bovine,Porcine
Conjugation	HRP
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Thrombin light chain: 328-363/622
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	4/64kDa
Uniprot ID	<b>P00734</b>
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.