

Product Datasheet

ITIH1 antibody {HRP} orb101727-HRP

Description	Rabbit Polyclonal antibody against ITIH1 conjugated to HRP.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ITIH1: 195-290/911
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	101kDa
Uniprot ID	P19827
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.