

Product Datasheet

RING2 antibody {CF594} orb1895-CF594

Description	Rabbit polyclonal antibody against RING2 conjugated to CF594.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Canine, Equine, Bovine, Porcine, Gallus
Conjugation	CF594
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human RING2: 51-160/336
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	100 µg in 500 µl (0.20 mg/ml)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	37kDa
Uniprot ID	Q99496 , Q4KLY4 , Q9CQJ4
Dilution Range	To be determined by end user
Expiration Date	12 months from date of receipt.