

## Product Datasheet

### Desmin antibody orb1151455

Description	Mouse monoclonal antibody to Desmin
Species/Host	Mouse
Reactivity	Bovine, Canine, Equine, Feline, Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Recombinant full-length human desmin protein was used as the immunogen for the Desmin antibody.
Preservatives	0.2 mg/ml in 1X PBS with 0.1 mg/ml rAlbumin (US sourced), 0.05% sodium azide
Storage	Aliquot the Desmin antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.
Note	For research use only.
Isotype	Mouse IgG2b, kappa
Clonality	Monoclonal
Purity	Protein A/G affinity
Formula	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide
Clone ID	DES/3255
Uniprot ID	<a href="#">P17661</a>
Dilution Range	Immunohistochemistry (FFPE): 2-4ug/ml
Expiration Date	12 months from date of receipt.
Application Notes	Optimal dilution of the Desmin antibody should be determined by the researcher.