


Product Datasheet

Tau-4 antibody orb1152101

Description	Tau-4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human Tau-F protein (1-441/441 aa)
Target	Tau-4
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	46/79 kDa
Uniprot ID	P10636
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane
1: Human
SH-SY5Y Cell; Lane
2: Human
MCF-7 Cell; Lane 3: Human
Jurkat
Cell; Primary: Rabbit
Anti-Tau-F
Protein Antibody at 1:500
dilution.