

Product Datasheet

Human Tau-4 antibody orb1152037

Description	Human Tau-4 antibody
Species/Host	Mouse
Predicted Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Tau-4 protein (1-441/441 aa)
Target	Human Tau-4
Form/Appearance	Liquid
Preservatives	0.01M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
MW	46/79 kDa
Clone ID	5A12
Uniprot ID	P10636
Dilution Range	WB=1:500-2000, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.