

Product Datasheet

Desmin antibody orb1151963

Description	Desmin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Desmin
Target	Desmin
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Recombinant
MW	52 kDa
Uniprot ID	P17661
Dilution Range	WB=1:500-2000, IHC-P=1:200-1:1000, ICC=1:50, IF=1:100-500
Expiration Date	12 months from date of receipt.