

## Product Datasheet

### Tick-borne encephalitis virus Genome polyprotein Antibody orb1152260

|                        |                                                                                   |
|------------------------|-----------------------------------------------------------------------------------|
| <b>Description</b>     | Tick-borne encephalitis virus Genome polyprotein Antibody                         |
| <b>Species/Host</b>    | Rabbit                                                                            |
| <b>Reactivity</b>      | Virus                                                                             |
| <b>Conjugation</b>     | Unconjugated                                                                      |
| <b>Form/Appearance</b> | Liquid                                                                            |
| <b>Concentration</b>   | 500 µg/ml                                                                         |
| <b>Storage</b>         | At -20°C for one year, at 4°C for one month. Avoid repeated freezing and thawing. |
| <b>Note</b>            | For research use only.                                                            |
| <b>Isotype</b>         | Rabbit IgG                                                                        |
| <b>Clonality</b>       | Polyclonal                                                                        |
| <b>Uniprot ID</b>      | <a href="#">Q01299</a>                                                            |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                   |