



Product Datasheet

Caspase 9 antibody orb10242

Description	Caspase 9 antibody
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Predicted Reactivity	Bovine
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Caspase-9 subunit p35 (11-120/416 aa)
Target	Caspase-9
Form/Appearance	Liquid
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only.
Isotype	IgG
Clonality	Polyclonal
MW	35/50 kDa
Uniprot ID	P55211
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ICC=1:100, Flow-Cyt=1?g/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin fixed and paraffin embedded rat brain tissue using Caspase 9 antibody