

## Product Datasheet

### ASK1/MAPKKK5 antibody orb11024

|                             |                                                                                  |
|-----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>          | ASK1/MAPKKK5 antibody                                                            |
| <b>Species/Host</b>         | Rabbit                                                                           |
| <b>Reactivity</b>           | Human                                                                            |
| <b>Predicted Reactivity</b> | Mouse, Rat                                                                       |
| <b>Conjugation</b>          | Unconjugated                                                                     |
| <b>Tested Applications</b>  | FC                                                                               |
| <b>Immunogen</b>            | KLH conjugated synthetic peptide derived from human ASK1 (321-460/1374 aa)       |
| <b>Target</b>               | ASK1                                                                             |
| <b>Form/Appearance</b>      | Liquid                                                                           |
| <b>Concentration</b>        | 1mg/ml                                                                           |
| <b>Preservatives</b>        | 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.            |
| <b>Storage</b>              | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| <b>Note</b>                 | For research use only.                                                           |
| <b>Isotype</b>              | IgG                                                                              |
| <b>Clonality</b>            | Polyclonal                                                                       |
| <b>MW</b>                   | 155 kDa                                                                          |
| <b>Uniprot ID</b>           | <a href="#">Q99683</a>                                                           |
| <b>Dilution Range</b>       | Flow-Cyt=3ug/Test, WB=1:500-2000, ELISA=1:5000-10000                             |
| <b>Expiration Date</b>      | 12 months from date of receipt.                                                  |

